



www.applied-maths.com

GeneMaths XT version 2.10 Upgrade - Release Note

31 May 2008

GeneMaths XT version 2.10 contains a number of bug fixes and technical improvements on the major upgrade version 2.0. Most improvements are therefore not directly visible as new features. Nonetheless it is highly recommended to install this upgrade as a more stable and mature version. For a full overview of new features in version 2.0, we refer to the 2.0 upgrade release note also available on the Applied Maths website.

Below is a list of new features in version 2.10:

- Additional algorithm to summarize value while using ***Collapse by group*** and ***Collapse by field***: in addition to **Average** now also **Weighted average** and **Median**.
- A ***Next array*** option in plot window to loop easily through different arrays while preprocessing
- Additional text field editing options: ***Add text field*** and ***Explode text field***
- Additional annotations ***Go Titles*** en ***KEGG titles***
- Annotations can optional be checked for newer verions on source web sites